

# ZBF150AC Blood Bank Refrigerator

## ELECTRIC - COMPRESSION WITH SURE CHILL® TECHNOLOGY




**ZERO**  
 A P P L I A N C E S

Application: Blood Bank Refrigerator  
 Type Of Appliance: Ice Lined (water lined) Refrigerator  
 Manufactured in: South Africa  
 Cold Chain Compliant: YES  
 IEC 60335 Compliant: YES

### Specifications

Gross Volume (L): 173  
 Drawers: 6  
 Blood Bag Storage: 72 x 450ml Blood bags  
 Hot Zone Appliance (43°C): Yes  
 Refrigerant: R134a (CFC Free)  
 Min Rated Ambient Temp: +5°C  
 Energy Source: Electric Mains 220-240V 50/60Hz or Electric Mains 115V / 60Hz Optional  
 Fuel And Cycle Type: Electric – Compression  
 Temperature Indicator: External Solar Thermometer



### Power Consumption

Power Requirements (as per ZLF150AC): 220 – 240V 50/60Hz. Included Is An ASP To Enable 150-280V  
 Power Consumption 43°C (as per ZLF150AC): 4 Hrs Mains Supply / 24 Hours  
 Power Consumption 32°C (as per ZLF150AC): 2½ Hrs Mains Supply / 24 Hours

### Performance Information

Hold Over Time: >110 Hours  
 Cooldown Time 32°C: Full Load Half Load  
 (Bloodbags Stable At 37°C): 19 Hours 11 Hours



\*Unit Uses **Sure Chill**® Technology. Water lining guarantees no freezing of stored items

### Accessories Optional Extras

Liquid PCM to mix with water  
 Drawers, 115V/60Hz nominal power option,  
 Plywood packaging, 2 padlocks

### Dimensions

Unboxed Dimensions (mm): (H) 1880mm x (W) 800mm x (D) 700mm  
 Shipping Volume (Crated): 1.22m³  
 Shipping Weight (Crated): 180kg

### Temperature Control

Temperature accuracy is exceptionally high; temperature varies by less than 1 °C throughout the refrigerated compartment during compressor cycling and holdover.

Zero Appliances (Pty) Ltd • Johannesburg, South Africa

 010 207 1600 • Int  +27 10 207 1600

 [info@zeroappliances.co.za](mailto:info@zeroappliances.co.za) •  [www.zeroappliances.co.za](http://www.zeroappliances.co.za)

E&OE Errors and Omissions Excepted